

Fersentlastungsorthese

nach Dr. Settner / OTM Münch

Heel relief orthosis




designed by Dr. Settner / CPO Münch

DE Anleitung für Patienten

EN Instruction for patients



DE Anlegen
EN Application





DE Klassischer Therapieverlauf

8. – 12. Tag		Volle Belastung ohne Gehstützen
4. Woche		Erste Röntgenkontrolle
6. Woche		Das erste Druckaufbaupolster wird eingesetzt
8. Woche		Zweite Röntgenkontrolle
		Das Zweite Druckaufbaupolster wird eingesetzt
10. Woche		Das dritte Druckaufbaupolster wird eingesetzt
11. Woche		Falls erforderlich, Maßnahmen für orthopädische Schuhe oder Einlagen
12. Woche		Abschluss der Behandlung

Der individuelle Therapieplan wird vom Arzt festgelegt!

EN Classical Therapy plan

Day 8 – 12		Full load without crutches
Week 4		First radiological check
Week 6		First compression pad is inserted
Week 8		Second radiological check
		Second compression pad is inserted
Week 10		Third compression pad is inserted
Week 11		If necessary, taking measurements for orthopaedic shoes or insoles
Week 12		Conclusion of treatment

Individual therapy is determined by the attending physician!



DE Arzt

EN Physician



DE Sanitätshaus

EN Specialist dealer

DE Reinigung der Orthese

EN Cleaning of the orthosis

Polster/
Pads



Orthese/
Orthosis



Otto Bock HealthCare GmbH
Max-Näder-Str. 15 · 37115 Duderstadt/Germany
T +49 5527 848-1706 · F +49 5527 72330
export@ottobock.de · www.ottobock.com