



ottobock.

Silicone cover for lower limb prostheses

Modern and custom-made

Quality for life



A perfect match – Custom-made silicone covers

Custom-made silicone covers for lower limb prostheses perfectly restore your external appearance.

The silicone cover, covers the prosthesis up to the knee joint and forms a harmonious part of your body image. The custom-made silicone cover is modelled to resemble your existing leg.

Thanks to the use of silicone, it is not only possible to recreate the colour of your skin surface perfectly, but also to capture the smallest details such as small veins, skin pigments or leg hair. This includes a custom-made reproduction of your toenails. You can also paint acrylic toenails.

If you opt for the "Natural" design, you can influence, prior to the final cosmetic design, the look of your silicone cover in order to obtain the best possible aesthetic appearance. Even special requirements such as tattoos are available.

Advantages of silicone:

- Original colour based on your existing leg
- Natural shape
- Robust material (water-repellent, UV-proof, temperature-resistant)
- Easy cleaning with water and soap
- Reduced sensitivity to colour transferring from clothes, e.g. jeans

Each silicone cover is unique, fabricated carefully by hand according to your requirements.

