

Liner Auswahlhilfe

Eigenschaften	Modelle													
Produkt	Skeo	Skeo Pure	Skeo Skinguard	Skeo 3D	Skeo Unique	Caleo	Caleo 3D	Uneo	Uneo Pure	Uneo Skinguard	Uneo Flex	Uneo 3D	Uneo Fresh	Uneo Unique
Oberschenkel			●	●	●									●
Unterschenkel	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Stumpf mit viel Weichteilgewebe	●	●	●	●	●									
Stumpf mit trockener Haut						●	●							
Knöcherner, narbiger Stumpf								●	●	●	●	●	●	●
Antibakteriell			●		(●)*					●		(●)*		(●)*
Mit Duftstoff													●	(●)*
Anatomisch geformt				●	●		●				●	●	●	●
Maßanfertigung					●									●
Transparent		●							●	(●)*	(●)*	(●)*	(●)*	(●)*

* optional