

## Ventus – unique in its diversity

---

### A reliable companion in all life situations

With the Ventus, Ottobock has successfully developed a wheelchair that meets the differentiated needs of active wheelchair users. It combines outstanding quality, ergonomics and functionality with an appealing design at a high level. At the same time, the Ventus meets the demand for mobility in all its facets at a very attractive price.

Quality of life means travel, sports or getting from A to B without help. But ordinary activities like riding in a car can quickly become an obstacle for wheelchair users. The Ventus is designed especially to meet such everyday challenges. Without impairing the driving stability of the wheelchair, the user can remove the wheels at any time in just a few steps and fold the backrest down onto the seat bottom. This makes the Ventus extremely compact for stowing. With a frame that weighs a mere eight kilograms, it is also a real lightweight. These user-friendly product characteristics make it possible for the driver to stow the wheelchair easily and independently, which means this reliable companion is always by the user's side.

#### **All advantages at a glance**

Based on many years of empirical values accumulated in the fabrication of high-quality mobility solutions, the Ventus is made from components that have already proven themselves in the market. The result is a permanently welded rigid frame wheelchair made of lightweight aluminium material, which is both robust and sturdy. As a result, it features outstanding driving stability and manoeuvrability along with efficient force transmission. The load capacity varies between 100 and 140 kilograms depending on the version. In addition to the higher load capacity, the Heavy Duty model features a seat width of 50 centimetres. Especially for users who are still growing, the Ventus is available with seat widths of 28-36 centimetres as well.

It is also suitable for initial fittings thanks to the height-adjustable back. The raised back offers greater support at the outset, which can be reduced step by step depending on the user's muscle development. An integrated shock absorber that dampens the rear axle helps prevent spasms and is suitable for users with high pain



Quality for life

sensitivity as well. The Ventus proves that it meets a wide variety of requirement profiles, and not only due to its ergonomic potential. Its appealing design is another of its strengths. The combination of a sporty frame structure with a palette of 14 trendy colours along with four modern Pure Metal hues makes it a true eye-catcher. From ice blue to candy red to deep black matte, something is offered to suit every taste.

**For further information:**

Karsten Ley

Director Corporate Communication

Otto Bock HealthCare GmbH, Max-Näder-Str. 15, 37115 Duderstadt

Phone: +49 (0)5527 848-3036, Fax: +49 (0)5527 848-3360

E-mail: [press@ottobock.com](mailto:press@ottobock.com)

Internet: [www.ottobock.com](http://www.ottobock.com)